

April 7, 2017



Product Safety Alert

## **SCIEX Recalls Varian TV 801 Turbo Pumps on the API 4000™, 4000 QTRAP®, and API 5000™ Mass Spectrometers Due to Mechanical Hazard**

The Electrical Safety Authority is notifying the public that SCIEX has announced a voluntary recall of the following commercial product. Users should stop using recalled products immediately and follow the steps indicated in the Recall notice.

Name of Product: API 4000™ LC/MS/MS, 4000 QTRAP® LC/MS/MS, API 5000™ LC/MS/MS For Research Use Only. Not for use in diagnostic procedures

Units: 3500 Manufacturer/Importer: MDS SCIEX, MDS Analytical Technologies AB SCIEX

Manufactured in: Canada, Singapore Hazard: Injury or property damage due to expelled parts

Incidents/Injuries: Two incidents have been reported where property damage was sustained due to expelled parts. No injuries reported.

Sold by: AB SCIEX LP 71 Four Valley Drive Vaughan, ON L4K 4V8 and its worldwide affiliates under a number of different brand names, including SCIEX, AB SCIEX, MDS SCIEX, and Applied Biosystems.

When distributed: 2002 to 2017

Remedy: Shut down the instrument, disconnect from the mains supply to prevent further use until the instrument is inspected and serviced by AB SCIEX.

Consumer Contact: AB SCIEX Canada and USA: Telephone: 1-877-740-2129 (toll-free) Other regions go to <https://sciex.com/about-us/contact-us> to find contact information. Visit the web site below for further information: <https://SCIEX.com/tv801safetynotice>.

Remarks: The recalled units are certified to Canadian Standards by CSA, certification number LR 81426. For more information about CSA product certification process please visit: [www.csagroup.org](http://www.csagroup.org).